
















PARTS LIST

ITEM		DESCRIPTION	PART NO.	QTY.	
1		BULKHEAD CONNECTOR	B-5568-000	1	
2		BULKHEAD COUPLER, VALVED	PA-1633-000	1	
3		UNVALVED NIPPLE	PA-1635-000	1	
4		HEX NUT, NYLON, 3/4-10	PA-1636-000	1	
5		COUPLER NUT	PA-1634-000	1	
6		WASHER, SEALING, NEOPRENE	PA-1637-000	2	
7		WASHER, FLAT, NYLON	PA-1638-000	2	
8		ELBOW, 22MM FEM X 7MM MALE	B-5603-000	2	
9		CAPLUG (12), WHITE	62905510W	1	For use with Silhouette only
10		CONNECTOR, 22MM M/F, 90 DEG ELBOW	1571-22	1	
11		HOSE, CLEAR W/WHITE BRAID, OXYGEN, 5/16 ID, 25FT LONG	PA-1640-001	1	
12		TUBING, 9/16 OD, 3/8 ID, PVC, CLEAR, 25 FT LONG	PA-679-CLR25	1	
13		TUBING, .625OD-3/8ID, GREY, 2" LONG	A-4484-000	1	
14		HOLE PLUG, 3/4", WHITE	PA-1639-000	2	Not shown in instructions
15		TIE WRAP 5.5" NYLON, WHITE	31-00156-00	6	Use as needed




Parker Hannifin Corporation
Porter Instrument Division
245 Township Line Rd.
Hatfield, PA 19440
(215) 723-4000
www.porterinstrument.com/dental

REMOTE CHAIR MOUNT KIT INSTRUCTIONS

PART NUMBER - 2025

COMPONENTS INSIDE CABINET



BELOW FLOOR

CAUTION: When installing hoses for the 2025 Chair Mount Kit ensure there are no kinks (or potential for kinking) as this will affect the performance of the system.

